

Red Excitable Fluorescent Dyes, Excitation, Emission & Uses

LASER

EX/EM nm

FLUOROCHROME

USE

Red	•633/660	Allophycocyanin, APC	antibody fluorochrome
Red	•633/672	AlexaFluor647	antibody fluorochrome
Red	•633/720	Alexa Fluor 700	antibody fluorochrome
Red	•633/720	APC-Cy5.5	antibody fluorochrome
Red	•633/780	APC-Cy7	antibody fluorochrome
Red	•633/780	APC-H7	antibody fluorochrome
Red	•633/660	Cy5	antibody fluorochrome

Red	•633/680	DRAQ5	live cell cycle analysis
Red	•633/780	DRAQ7	cell viability
Red	•633/680-730	Vybrant DyeCycle Ruby	live cell cycle
Red	•633/660	FxCycle Far Red	fixed cell cycle
Red	•633/660	To-Pro3	cell viability, fixed cell cycle
Red	•633/660	DiIC(5)	mitochondrial membrane potential
Red	•633/660	LysoTracker Deep Red	lysosomal mass

Red	•633/740 •633/775	Zombie NIR APC-eFluor780	Fixed cell viability antibody fluorochrome
Red			
Red	.		